

# Angelman syndrome (46,XX,del(15)(q11q13)) Reference Standard

### **CBPJ0006**

## I. Description

NIPT reference products use gDNA from human cell lines as raw materials, which are fragmented and then added to human plasma to simulate real clinical samples. They can monitor the quality control of sample extraction, library preparation and sequencing of non-invasive prenatal testing (NIPT) kits. They are suitable for the detection of cell-free DNA (cfDNA) in peripheral blood of pregnant women using the whole genome method (massive-parallel sequencing MPS).

CB-Gene NIPT reference products can be customized for common chromosomal abnormalities: aneuploidy of four common chromosomes (9/21/13/18); aneuploidy of one sex chromosomes is aneuploid (XXY/47); Common microdeletion and microduplication syndrome types, such as DGS, AS, PWS, 11q23.3 del, etc.

#### II. General information

Name	Angelman syndrome	
	(46,XX,del(15)(q11q13)) Reference	
	Standard	
Cat. No.	CBPJ0006	
Format	Genomic DNA	
Unit Size	1ug	
Intended Use	Research Use Only	
Buffer	Tris-EDTA	
Storage Conditions	2-8°C	
Expiry	36 months from the date of manufacture	



# **III.**Technical Data

Mutation Angelman syndrome

Karyotype 46,XX,del(15)(q11q13)

Chromosome

(GRCh37)

abnormal region

chr15:22,770,422-28,644,578

(CMA)

Mutation Type Loss

Length (bp) 5874157

ACMG variant


Pathogenicity classification

# **IV.Representative Data**

#### Z-Score:

chr	Z-Score	chr	Z-Score
chr1	4.56	chr12	1.13
chr2	2.78	chr13	-1.84
chr3	-0.21	chr14	0.01
chr4	-3.05	chr15	-8.52
chr5	3.15	chr16	0.96
chr6	-1.16	chr17	1.51
chr7	3.19	chr18	-4.28
chr8	-2.9	chr19	3.31
chr9	-0.23	chr20	-0.88
chr10	-2.23	chr21	-2.38
chr11	-1.14	chr22	0.35









## Chromosome simulated karyotype:

